



- SMP-PAD IS MADE FROM SMP(SHAPE MEMORY PLASTIC).
- SMP-PAD STARTS TO CHANGE THE SHAPE AT BODY TEMPERATURE (36 DEGREES).
- 3 SIZE, S/M/L AVAILABLE
- AVAILABLE WITH ONLY BOX TYPE(SCREW TYPE)
- SMP-PAD WILL BE SUITABLE FOR COMFORTABLE FITTING OF NOSE PADS.
- SMP(SHAPE MOMORY PLASTIC) IS USED IN MEDICAL AND HEALTH CARE FIELD.

NOTE: STORE SMP IN COOL TEMPERATURE TO AVOID ACTIVATION.

SMP STARTS TO CHANGE AT BODY TEMPERATURE.





